

## AuDX Pro The Preferred Portable Diagnostic OAE System for All Ages



### Flexible System for all Practice Types

- The AuDX Pro is the ideal screener for use in UNHS programs, primary care physician's offices and schools
- The AuDX Pro II is ideal for Audiologists and ENT Physicians or anyone who needs increased system flexibility
- The AuDX Pro Plus combines portable screening and full computer based diagnostic OAE capabilities for ultimate flexibility

*All infants should have their hearing screened before 1 mo of age. Diagnosis of hearing loss should be completed by 3 months of age so that appropriate rehabilitation can begin before the infant is six months old.*

### Why professionals choose AuDX Pro:

- Objective test – no patient response required
- Flexible testing options with DPOAE or TEOAE
- Pass/Refer result provided – no interpretation required
- Easy to use with simple workflows
- Fast, accurate results – as fast as 10 seconds per ear
- Preprogrammed and customizable protocols
- Adult screening protocol – avoids over-referral issues associated with pediatric protocols
- 12 frequency diagnostic testing ability
- Large color display – enhanced operation & data viewing
- Label printer – results documented in seconds
- Operates on rechargeable battery or AC power – no down time during recharge
- Backed by over 30 years of experience with hearing screening and diagnostic product design and manufacturing
- 24/7 Technical Support

# The Portable Diagnostic OAE Hearing System of Choice for Providers of World Class Care

## AuDX Pro

- Standard system includes
  - 50 test memory
  - **DPOAE or TEOAE** screening protocols
- Options include:
  - Scout Link software for data transfer to PC
  - HATS data management software
  - Additional (3) DP protocols including "adult protocol"
  - Protocol Setup Software
  - Label printer
  - Large carry case (holds printer and AuDX Pro)

## AuDX Pro II

- Standard system includes:
  - 100 test memory
  - **DPOAE and TEOAE** screening protocols
  - Additional (3) DP protocols including "adult protocol"
  - Protocol Setup Software
  - 12-frequency testing within single protocol
- Options include:
  - Scout Link software for data transfer to PC
  - HATS data management software
  - Label printer
  - Large carry case

## AuDX Pro Plus

- Standard system includes:
  - 100 test memory
  - **DPOAE and TEOAE** screening protocols
  - Additional (3) DP protocols including "adult protocol"
  - Protocol Setup Software
  - 12-frequency testing within single protocol
  - Scout Collection software (requires PC) for full diagnostic OAE
  - Scout Link software for data transfer to PC
- Options include:
  - Label printer
  - Large carry case

### HATS™


Allows download of AuDX data into this powerful, PC-based screening data management software. Generate patient reports, statistical reports for quality assessment and export data into a format compatible with other 3rd party databases (requires a PC).

### Scout Link

Allows download of AuDX data into a reader version of the Scout PC-based application for reviewing, storing and printing data on the PC (requires a PC).

Easy to clean,  
durable Ear Probe

Your Local Distributor:



*Bio-logic supplies provide quality & value while ensuring optimal performance*



Premie Tree Tip®

Tree Tip®

Pediatric Tip

Pediatric Foam Tip

Adult Foam Tip

Jumbo Foam Tip

**natus**

Natus Medical Incorporated  
www.natus.com

Visit our NERVE Center® education portal at [nervecenter.natus.com](http://nervecenter.natus.com)

© 2015 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.