



Natus® SleepWorks™

Comprehensive, flexible and secure sleep diagnostic software

Natus SleepWorks is a comprehensive, flexible and cyber-secure sleep diagnostic solution that simplifies the setup, collection, review, analysis and reporting of comprehensive sleep studies. With advanced HIS/EMR integration and unmatched reliability, SleepWorks is a complete diagnostic package designed to meet the needs of private practice clinics, sleep clinics, hospitals and academic medical centers.

Enhanced workflow for sleep studies

- Ongoing study parameters at the forefront of view in a configurable Heads Up Display
- Configurable observation notes – automate tech notes with custom form fields and flexible reporting
- Analyze and review an ongoing studies, including video
- Capture PAP treatment and supplemental O₂ levels easily using customizable therapy types that simplify titration and reporting

Intuitive Clinical Experience

- Powerful event detection, trending and analysis options with fully synchronized HD video
- Record data directly to a server and access studies from any networked station
- Reduce size of archived files using advanced clipping and pruning features
- Use the Natus Database to quickly locate patient studies based on a variety of custom parameters



What's new in SleepWorks

- ★ Windows® 10 Support
- ★ New hardware platform
- ★ Improved central administrator functions
- ★ Automated observation notes
- ★ EDF/EDF+ import and export

SleepWorks provides unmatched reliability and fault tolerant distributed networking solutions

Robust Security Features

- Meets HIPAA requirement via additional XLSecurity layer
- Define access to features, studies and programs based on logged-in user or role
- Automated power recovery
- Microsoft® BitLocker® drive encryption support

Site Management tools for enhanced productivity

- Enhanced security features with Active Directory integration
- Credentials cache system
- Automated site-wide software deployment with optional Cerebrum Enterprise Solution

A new cutting-edge hardware platform bringing together long standing reliability and advanced technology

Natus' new hardware platform is inspired by a rich heritage of expertise, drawing from our top sleep brands. The amplifiers consist of two main parts, a breakout box and a wall or cart mounted universal Natus Base Unit.

Breakout Box Features

- Small, rugged and lightweight amplifiers
- Powered sensor inputs
- Integrated pressure sensor
- Integrated pulse oximeter
- Up to 103 channels

Natus Base Unit Features

- Touchscreen display for real time interaction with ongoing studies
- Visual impedance checking and adding notes
- Perform bio-calibrations while in the patient room
- Universally compatible with several Natus EEG/PSG amplifier breakouts

Full range of powerful amplifiers:



Standard sleep testing with **Embla SDx**



Advanced neuro-sleep monitoring with **Embla NDx**



Now integrated with **Embletta MPR Scalable Sleep System**

natus
neuro

Natus Neuro

3150 Pleasant View Road
Middleton, WI 53562 USA
Tel: 1-800-356-0007
1-608-829-8500
Fax: 1-608-829-8709

www.natus.com

To learn more about Natus Neuro products contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**